



SAFETY SYSTEM®
YOUR PROTECTOR

Sir Safety System S.p.A. Single-member company
Wholly-owned subsidiary of Sir Holding Srl
Via dei Fornaciaci, 9 - Santa Maria degli Angeli
06081 Assisi - Perugia - Italy
Tel.+39 075.8043737 - Fax +39 075.8048686

VAT Nr. IT 03359340548
Viale Europa 74
20090 Cusago Milano
Tel. +39 02.90394575
Fax +39 02.90394441

TECHNICAL SHEET

export@sirsafety.com
www.sirsafety.com



UNI EN ISO 9001:2015
UNI EN ISO 14001:2015
UNI EN ISO 45001:2018

Name: MEDIBAG 1 first aid case	Code: FH1213	FH12 MEDICAL-HEALTH EQUIPMENT - FIRST AID
---------------------------------------	---------------------	--

Description

First aid case in impact resistant polypropylene. Fitted with wall attachment plate and handle for easy carrying. Supplied with first aid items required by It. M.D. 388 15/7/03 attachment 2 and It. L.D. 81 9/4/08 art. 45. Suitable for premises up to two employees.
Dimensions: 25 x 19 x 9 cm.

Content:

- no. 1 copy M.D. 388 15/7/03
- no. 2 pairs of sterile gloves
- no. 1 bottle of antiseptic solution, 125 ml, 10% Povidone-iodine
- no. 1 sterile saline solution, 250 ml
- no. 1 pack of sterile compressed gauze 18x40 cm
- no. 3 packs of sterile compressed gauze 10x10 cm
- no. 1 sterile forceps
- no. 1 pack cotton wool
- no. 1 Plastosan box of 10 assorted plasters
- no. 1 adhesive plaster roll 5 m x 2.5 cm
- no. 1 gauze bandage 3.5 m x 10 cm
- no. 1 pair of Lister scissors 14.5 cm
- no. 1 tourniquet
- no. 1 instant Ice Pack, disposable
- no. 1 sanitary waste bag 250x350 mm

Supplied with first aid items required by It. M.D. 388 15/7/03 attachment 2 and It. L.D. 81 9/4/08 art. 45.



Category	Colours	Available Sizes In stock: U
Quantity/Carton 10	EC Standards	Important Information Suitable for working environments with up to two employees. Read carefully the user sheet included. Make sure the item is in good condition, undamaged before use. Once you become aware of a damage, you should immediately replace it. Keep in a dry place, far from heat sources.
Dexterity		Storage and cleaning