

Sir Safety System S.p.A. Single-member company Wholly-owned subsidiary of Sir Holding Srl Via dei Fornaciai, 9 - Santa Maria degli Angeli 06081 Assisi - Perugia - Italy Tel.+39 075.8043737 - Fax +39 075.8048686 VAT Nr. IT 03359340548 Viale Europa 74 20090 Cusago Milano Tel. +39 02.90394575 Fax +39 02.90394441

TECHNICAL SHEET

export@sirsafety.com www.sirsafety.com





Name: Ladies' lab coat Code: MC2646K1 **MC26 PROFESSIONAL CLOTHING** Description Ladies' lab coat made of 65% polyester, 32% sanforised cotton and 3% lycra, 150 g/m, featuring a womencs fit, knee-length, two lower pockets and chest pocket with pen holder, concealed six-button closure, cuffs fastened with buttons and rear martingale. - Suitable in warm environments. - Easy to wash and sanitise. Protection against risks of dirt and detrimental actions of slight damage. The product is designed and manufactured to comply with the requirements of EU Regulation 2016/425 and subsequent amendments. V@-Á, æn∄ Á¢¢o^¦}ænÁæà¦a‰Áa ÁJÒSUË/ÒÝí ÁÙæa)åælåÁF€€Á OEKO TEX® certified Available Sizes Category Colours I WHITE In stock: 38 - 56 Quantity/Carton EC Standards Important Information EN ISO 13688 Suitable in medical laboratories and as visitor's coat. Store in a cool, dry place, away from heat sources and away from light. 30 Before each use, carry out a visual check to ensure that the device is in perfect condition. undamaged and clean, and replace it ifit is not intact. Follow the recommended washing instructions. Dexterity Storage and cleaning ¥@7 🖄 🕅 🕑