

Name: Men's lab coat

Sir Safety System S.p.A. Single-member company Wholly-owned subsidiary of Sir Holding Srl Via dei Fornaciai, 9 - Santa Maria degli Angeli 06081 Assisi - Perugia - Italy Tel.+39 075.8043737 - Fax +39 075.8048686 VAT Nr. IT 03359340548 Viale Europa 74 20090 Cusago Milano Tel. +39 02.90394575 Fax +39 02.90394441

Code: MC2647K1

TECHNICAL SHEET

export@sirsafety.com www.sirsafety.com





Description
Men's lab coat made of 65% polyester and 35% sanforised cotton, 150 g/m . Knee-length, two lower pockets and chest pocket with pen holder, concealed five-button closure, cuffs fastened with buttons and rear martingale.
Suitable in warm environments.
Easyto wash and sanitise.
Protection against risks of dirt and detrimental actions of slight damage.

The product is designed and manufactured to comply with the requirements

of EU Regulation 2016/425 and subsequent modifications.

V@A(, æn), Ár¢c^\}æµÁæaà¦a3a/≨rÁUÒSUË/ÒÝí ÁÙæa) 忦å/aF€€Á certified



MC26 PROFESSIONAL CLOTHING

| certified | | IEX® |
|-----------------|--------------|---|
| Category | Colours | Available Sizes |
| I | WHITE | In stock: 44 - 62 |
| Quantity/Carton | EC Standards | Important Information |
| | EN ISO 13688 | Suitable in medical laboratories and as visitor's coat. |
| 30 | | Store in a cool, dry place, away from heat sources and away from light. Before each use, carry out a visual check to ensure that the device is in perfect condition, undamaged and clean, and replace it ifit is not intact. |
| | | Follow the recommended washing instructions. |
| Dexterity | | Storage and cleaning |
| | | |

ΟΕΚΟ